



LEGRAND BED EXTENSION - TO BE ASSOCIATED CU CAT.NO 0 766 06 OR 0 766 07

Caracteristici tehnice

Antimicrobial door units for rooms

For installation at the entrance to the room, can be seen from the bed

Connect directly to the BUS/SCS cable Cat.Nos 0 492 32/33

With indicator and pushbutton for calls, nurse presence and 2nd nurse presence

Comply with standards VDE 0834-1 and 0834-2

Fault display

4-module wiring mechanisms supplied complete

Flush mounting in 50 mm deep box or surface mounting with frame Cat.No 0 766 14

Can be acknowledged by a pushbutton with indicator Cat.No. 0 782 04

Antimicrobial : contains silver ions which prevent the development of surface bacteria

Bed extension

To be associated with Cat.No 0 766 06 or 0 766 07

Allows the connection of 4 components to 1 call unit DIN, 6 module

Programmation via virtual configuration

Caracteristici generale

BUS/SCS nurse call unit - signalling and acknowledgement units for rooms

BUS/SCS system with star wiring around the door unit:

call centralisation

acknowledgement

communication with nursing control unit and all elements of the system via BUS/SCS

signalling forwarded to all indicators (e.g.: corridor overhead light unit, door unit, hand-held remote control unit)

Pret: 1.391,79 LEI (TVA inclus)

Detalii online: <https://www.materialelectrice.ro/bed-extension-to-be-associated-cu-cat-no-0-766-06-or-0-766-07-222903>